INSTRUMENTAL QUIRÚRGICO

FABRICADO POR /MANUFACTURED BY: IMECO S.A.

Roca 2054 José León Suárez

Provincia de Buenos Aires

☎ (5411)4729-1110

✉: info@imeco.com.ar

www.imeco.com.ar

Industria Argentina

Made in Argentina

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Reutilizable

Conserve a temperatura ambiente

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Read carefully the User's Manual.

Description:

IMECO S.A Instrumentation consists of manual surgical instruments for use during hip, knee, shoulder, elbow or wrist arthroplasty, spinal fusion or trauma surgeries.

Utilization:

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.

- IMECO S.A trial components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus preserving the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- The IMECO S.A Surgical Protocols provide additional procedural information.

Warnings:

- Care should be taken not to cut through surgical gloves when handling any sharp-edge orthopedic device.
- Due to different manufacturer's employ differing design parameters, varying tolerances, different materials and manufacturing specifications, IMECO S.A Instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of IMECO S.A for the performance of the resulting implant.
- Malfunctions of instruments, including, but not limited to, weld fracture, thread gauling and material fatigue have occurred in a small percentage of cases.

Adverse effects:

- Serious complications may be associate with any joint replacement, spinal or trauma surgery. These complications include, but are not limited to: genitourinary disorders, vascular disorders, including thrombus, bronchopulmonary disorders, including emboli, myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotrophic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of surgical instruments and the impaction of the component into the preparation.

Instrumental Care:

- Before any kind of surgery, the material to be used must undergo sterilization and a cleaning process.
- During surgical and prosthetic procedures, the instrumental must be used only for what it is meant to. Any other usage of the instruments may cause fractures, deformations and/or diminish the lifespan of these.
- Immediately after surgery, any blood, bone or tissue residue must be removed. Do not leave on any surface. Always dismount compound instruments in order to guarantee the removal of all possible residues.
- As quickly as possible, proceed to disinfect and clean the surgical material by using adequate products, carefully watching the manufacturer's instructions in relation to the dose, action time and temperature.

Cleaning must be performed after disinfection, preferably using ultrasound equipment.

Manual Cleaning:

Its main objective is the removal of organic matter from the instrumental. This procedure must be initiated as soon as possible.

- Only use nylon brushes. Never use steel sponges or brushes.

- Use demineralized or distilled water with neutral soap or enzymatic detergent.
- Rinse strongly in common water and immediately dry.

Ultrasound Cleaning:

- Place the instrumental in adequate amounts to prevent damage.
- Never throw the instruments in the ultrasound cleaning bucket. These must be carefully placed inside the bucket to prevent material damage.
- Always remember to remove any cleaning product residues in common water and dry clean the instrumental immediately thereafter.

Important Notice:

- Aluminum pieces must not be placed in ultrasound equipment.
- Never use products with vast quantities of bleach, oxidizing substances, acids or any salts in general (i.e. hypochlorite sodium substances, oxygenated water, physiological serums, etc.). If any of these products are used, this may incur in localized or contact corrosion.
- Any association of pressurized, moist, high-temperature chemical or biological residues will cause corrosion/oxidation of the instruments.

Dry Cleaning / Inspection:

After the cleanup process, all instruments must be treated with ethylic alcohol and left to dry. Instruments with articulations or cavities that may difficult the drying process must be dried with a blast of compressed air (controlled process). Instruments must not be sterilized while damp. Prior to sterilizing, a cautious inspection of the instrumentals must be conducted in order to make sure there are no possible damages (corrosion, loss of edge, etc.), as well as any organic matter residues that were not removed.

Special Cares:

- Never let the instrumental to air dry. The permanence of moisture on metal could cause stains and corrosion due to the presence of minerals and several other elements in water.

Sterilization:

The sterilization of surgical instruments must always be performed in any of the two following methods:

Method: Water Steam (Temperature: 121 ° C - Exposure time: 30 minutes)

Method: Ethylene Oxide - ISO 1135/1994 IRAM 37008

Prior to sterilization, the instrumental (already disinfected, clean and dry) must be properly conditioned.

Important Notice:

- Corroded or rusty instrumental pollutes the circuit inside the autoclave with particles, which can affect integral instrumental, causing stains and deterioration. For this reason,

always proceed in the cleaning and maintenance of the autoclave according to the manufacturer's specifications.

- Water, even when drinkable, has high concentrations of salts and chemical elements which deteriorate the instruments. Water used in the autoclave should be distilled, demineralized or deionized to avoid unwanted stains on said instruments.

DIT.IDIT.67

REV. 0

MAY/14